

1. Presentation

SILICON TUBING is a medical device used to connect a collection jar to a suction device or to connect a collection jar to a suction connector (probe, suction control valve, etc.). It can be connected to Ø9 mm nipples.

2. Versions



SILICON TUBING is available in the following versions:

- inside Ø = 7 mm x outside Ø = 13 mm sold in 25 m rolls
- inside Ø = 8 mm x outside Ø = 14 mm sold in 25 m rolls

3. Technical data

- Complies with standard ISO 10079-3: 2014.
- Lot number is specified on the packaging label.
- Ambient temperatures for use: -60°C to 200°C.

4. Symbols

	Manufacturer	 YYYY-MM	Use until YYYY (year) - MM (month)	LOT	Lot Nr
--	--------------	--	---------------------------------------	------------	-----------

5. Safety advice

- SILICONE TUBING should only be used by qualified medical personnel.
- Once cut, the length of SILICON TUBING must be at least 1.3 m long.
- Always check the use-by date specified on the packaging label.
- SILICON TUBING must never come into direct contact with the patient.
- SILICON TUBING is suitable for use in an MRI environment.
- To discard SILICON TUBING at the end of its service life, get in touch with an authorized organism. Do not dispose in bin for domestic waste.

6. Cleaning and disinfection

SILICON TUBING must be disinfected or sterilized after each patient. It can withstand 30 disinfection or sterilisation cycles.

SILICON TUBING can withstand autoclaving at 134°C for 18 min.

7. Transport and storage

- Store at between -10 and +40°C in a clean and dry place.
- Keep the packaging with the instructions until the product reaches the end of its service life.

8. Maintenance

No preventive or corrective maintenance is required.

In the event of damage or malfunction, return the device to the manufacturer or to an approved repairer.

9. Guarantee

- The device is covered by a 1-year warranty which does not cover damage or incidents caused by negligence, misuse, incorrect supervision or maintenance.
- Service life: 5 years (specified on the packaging label).
- Date when CE marking first affixed: 2000.